


pepette – The New Pipettor Generation

A lightweight precision instrument that withstands the rigors of daily laboratory life!



Superb Quality

The PEQpette – the perfect combination of highest quality materials

Durability and reliability

Pistons made of stainless steel or chemical-resistant polymer, sturdy teflon gaskets, plus shaft and casing made of chemical- and scratch-resistant material

Ruggedness

High-grade steel tip ejector, rustproof and deformation-resistant. Sturdy material, UV-resistant and autoclavable

Simple decontamination

Entire pipette is autoclavable and resistant against UV radiation. Pipette disassembly not required! Calibration is maintained even after repeated autoclaving.

Very Easy Handling

- Easy tip pick-up and ejection requiring only minimal effort
- Single-handed setting of volume via adjustment wheel or operating knob
- Flexibility through continuously adjustable volume
- Large, easy to read display
- Simple pipette calibration and maintenance via universal bung wrench
- Provides optimal fit for many of the most commonly available tips

Ergonomically Optimized

- Lightweight design and balanced weight distribution
- Consistent piston pressure prevents premature fatigue
- Can be used for both left- or right-handed individuals

peqlab 

Creating the future together.

Our Pledge

Our pledge is to provide excellent products at a fair price and to support our customers with a dedicated and experienced staff.



Ordering Information				
Item	Description	Accuracy (%)	Precision (%)	Order No.
PEQpette Single Channel Pipettes				
PEQpette 2E	Single Channel Pipette 0.1 – 2 µl	± 12 to ± 1.5	≤ 6.0 to ≤ 0.7	91-P3960-2
PEQpette 10E	Single Channel Pipette 0.5 – 10 µl	± 4.0 to ± 0.5	≤ 2.8 to ≤ 0.4	91-P3960-10
PEQpette 20E	Single Channel Pipette 2 – 20 µl	± 3.0 to ± 0.8	≤ 1.5 to ≤ 0.3	91-P3960-20
PEQpette 100E	Single Channel Pipette 10 – 100 µl	± 1.6 to ± 0.8	≤ 0.8 to ≤ 0.2	91-P3960-100
PEQpette 200E	Single Channel Pipette 20 – 200 µl	± 1.2 to ± 0.6	≤ 0.6 to ≤ 0.2	91-P3960-200
PEQpette 1000E	Single Channel Pipette 100 – 1000 µl	± 0.9 to ± 0.6	≤ 0.4 to ≤ 0.15	91-P3960-1000
PEQpette Multichannel Pipettes				
PEQpette 10M	8 Channel Pipette 1 – 10µl	± 8.0 to ± 2.0	≤ 6.0 to ≤ 1.2	91-P4608-10
PEQpette 50M	8 Channel Pipette 5 – 50 µl	± 4.0 to ± 1.6	≤ 2.5 to ≤ 0.6	91-P4608-50
PEQpette 200M	8 Channel Pipette 20 – 200 µl	± 3.0 to ± 1.0	≤ 1.5 to ≤ 0.6	91-P4608-200
PEQpette Carousel Stand				
Carousel stand	for 6 Single Channel PEQpettes			91-P3985

Further pipette models on request

